



**SUPERSONIC**  
imagine



# PRIOR COMMUNICATION


---

NEW TOOLS FOR LIVER ASSESSMENT

# SUPERSONIC IMAGINE – Pioneers in Innovative Imaging Modes

Pioneers in the development of **Innovative Imaging Modes** in response to clinical challenges, leveraging **unique UltraFast™ Imaging** performance.

... in 2009, SuperSonic Imagine introduced **ShearWave™ Elastography (SWE)** in the ultrasound market, as a unique solution for **tissue stiffness quantification**.



\* Imaging modes / features not yet available on SuperSonic Imagine products.

# NEW TOOLS FOR LIVER ASSESSMENT

Coming this Summer 2019 for your Aixplorer  
MACH 30 systems!

---

## NEW TOOLS FOR LIVER ASSESSMENT



---

Three new quantitative ultrasound tools, based on SSI-specific UltraFast™ technology, for the evaluation of different indicators of the chronic liver diseases severity (such as NAFLD, NASH, HBV, HCV...):

- **Ultrasound Beam Attenuation**
- **Speed of Sound**
- **Viscosity Imaging**




# NEW TOOLS FOR LIVER ASSESSMENT

## ❖ Attenuation


The attenuation measurement tool allows the quantification of the ultrasound beam attenuation in the liver through a B mode region of interest :

- Measurement of the decrease in amplitude of the ultrasound waves as they propagate, as a function of frequency.
- Measurement performed through a large ROI.
- Instantaneous result displayed on the 2D image (dB/cm/MHz), simultaneously acquired with the Speed of Sound thanks to UltraFast™ technology.
- The higher the attenuation coefficient; the higher the intrahepatic fat level ; the more severe grade of hepatic steatosis.\*

➔ **Attenuation tool seems to be an essential non-invasive biomarker to assess hepatic steatosis severity.**



Hepatic Steatosis – B Mode



Hepatic Steatosis with Attenuation quantification

# NEW TOOLS FOR LIVER ASSESSMENT




## ❖ Speed of Sound (SoS)


The Speed of Sound measurement tool allows the quantification of the intrahepatic speed of sound from B mode :

- SoS optimization in a targeted region of the liver, performed on a large amount of images (using UltraFast™ technology) to obtain the SoS giving the best image quality in this ROI. This SoS represents the speed of sound in the rest of the liver.
- Instantaneous result displayed on the 2D image (m/s), simultaneously acquired with Attenuation thanks to UltraFast™ technology.
- The lower the speed of sound ; the higher the intrahepatic fat level ; the more severe grade of hepatic steatosis.\*\*

→ SoS tool seems to be an essential non-invasive biomarker to assess hepatic steatosis severity.



Healthy Liver with Speed of Sound quantification



Hepatic Steatosis with Speed of Sound quantification

\*\*Ducquard B, et al. Ultrafast Adaptive Sound Speed Estimation for the Diagnosis and Quantification of Hepatic Steatosis: A Pilot Study. *Ultrasound Med Biol* 2018; No. 5. \*\*Imbaji et al. Robust sound speed estimation for ultrasound-based hepatic steatosis assessment. *2017 Phys. Med Biol*



# NEW TOOLS FOR LIVER ASSESSMENT



## ❖ Viscosity Imaging

This tool allows real-time visualization and quantification of tissue viscosity:

- Shear wave dispersion spectroscopy (analysis of shear wave propagation velocity at several frequencies), giving access to the shear wave dispersion measurement correlated to tissue viscosity (Pa.s).
- Combined in real-time with elasticity imaging (SWE), without compromising on 2D image quality.
- The higher the viscosity value; the higher the inflammation level (necro-inflammatory activity) in the liver.\*
- **Viscosity Imaging seems to be an essential non-invasive biomarker to assess chronic liver diseases severity.**



Healthy Liver



Courtesy of Pr V. Vilgrain – Beaulieu Hospital

Assumed Liver Cirrhosis


\*Christopher et al. Shear wave dispersion measures liver stiffness. *Ultrasound in Med & Biol.* 2011 \*Katafuchi Sugimoto et al. Viscoelasticity measurement in rat livers using shear wave elastography. *Ultrasound in Med & Biol.* 2008 \*Barnes et al. Shear Wave Dispersion: A Noninvasive Biomarker for Liver Fibrosis. *J Ultrasound Med* 2015

\*Delieux T et al. Shear Wave Spectroscopy in Vivo Quantification of Human Soft Tissue Viscoelasticity. *IEEE Transactions on Medical Imaging* 2016



# A COMPLETE SOLUTION FOR THE DIAGNOSIS OF CHRONIC LIVER DISEASES

- ❖ The arrival of these new innovative modes reinforces the positioning of SuperSonic Imagine as a major player in multiparametric ultrasound imaging of chronic liver diseases assessment, by obtaining the following information for the same patient:
  - Morphological characteristics through an amazing B-mode image quality.
  - Liver fibrosis staging through stiffness measurement (SWE™).
  - Hepatic Steatosis level through speed of sound measurement (SoS), and ultrasound beam attenuation measurement.\*
  - Necro-inflammatory activity level through viscosity measurement.\*
  - Perfusion (CEUS) and micro-vascularization (Angio PL.U.S) information allowing focal liver lesions characterization.
  - Information allowing prevention portal hypertension complications through SWE Spleen measurement.




\*Must be confirmed by clinical studies



# CLINICAL CASES

# CLINICAL CASE 1: HEALTHY PATIENT

Parameters	Results	Pathology Scoring
Attenuation	0.36 dB/cm/MHz	Probable S0 *
SoS	1595 m/s	Probable S0 **
Stiffness	3.47 kPa	F0/F1 ***
Viscosity	1.26 Pa.s	Probable A0 / S0 ****



Healthy Liver with Speed of Sound and Attenuation quantification



Healthy Liver with Stiffness and Viscosity visualization and quantification

\* Fujiwara et al., The b-mode image-guided ultrasound attenuation parameter accurately detects hepatic steatosis in chronic liver disease, *Ultrasound Med. & Biol.*, Vol. 00, No. 00, pp. 1-10, 2018


\*\* Dioguardi Burgio et al., Ultrasonic Adaptive Sound Speed Estimation for the Diagnosis and Quantification of Hepatic Steatosis: A Pilot Study, *Ultraschall Med.* 2018 Nov 5.

\*\*\* MKG.EC.100 (Rev B) - Aixplorer Aixplorer Ultimate\_Guide assessment fibrosis steatosis - Proposed diagnostic cut-off values and their performance for the evaluation of liver fibrosis severity by SWE™ - NASH - SuperSonic Imagine

\*\*\*\* Preliminary Internal R&D Study "Etude RND SSI hépatopathies chroniques 29072018" - SuperSonic Imagine

# CLINICAL CASE 2: OBESE PATIENT WITH KNOWN TYPE II DIABETES

Parameters	Results	Pathology Scoring
Attenuation	0.83 dB/cm/MHz	Probable S3 *
SoS	1495 m/s	Probable S1-S3 **
Stiffness	4.79 kPa	F0/F1***
Viscosity	1.88 Pa.s	Probable A1 to A2 / S2****



\* Fujiwara et al., The b-mode image-guided ultrasound attenuation parameter accurately detects hepatic steatosis in chronic liver disease, *Ultrasound Med. & Biol.*, Vol. 00, No. 00, pp. 1-10, 2018



\*\* Dioguardi Burgio et al., Ultrasonic Adaptive Sound Speed Estimation for the Diagnosis and Quantification of Hepatic Steatosis: A Pilot Study, *Ultraschall Med.* 2018 Nov 5.

\*\*\* MKG.EC.100 (Rev B) - Aixplorer Aixplorer Ultimate\_Guide assessment fibrosis steatosis - Proposed diagnostic cut-off values and their performance for the evaluation of liver fibrosis severity by SWE™ - NASH - SuperSonic Imagine

\*\*\*\* Preliminary Internal R&D Study "Etude RND SSI hépatopathies chroniques 29072018" - SuperSonic Imagine

# CLINICAL CASE 3: OBESE CHILDREN WITH KNOWN TYPE II DIABETES

Parameters	Results	Pathology Scoring
Attenuation	0.72 dB/cm/MHz	Probable S3 *
SoS	1450 m/s	Probable S1-S3 **
Stiffness	5.70 kPa	F0/F1 ***
Viscosity	1.92 Pa.s	Probable A1 to A2 / ****



\* Fujiwara et al., The b-mode image-guided ultrasound attenuation parameter accurately detects hepatic steatosis in chronic liver disease, *Ultrasound Med. & Biol.*, Vol. 00, No. 00, pp. 1-10, 2018



\*\* Dioguardi Burgio et al., Ultrasonic Adaptive Sound Speed Estimation for the Diagnosis and Quantification of Hepatic Steatosis: A Pilot Study, *Ultraschall Med.* 2018 Nov 5.

\*\*\* MKG.EC.100 (Rev B) - Aixplorer Ultimate\_Guide assessment fibrosis steatosis - Proposed diagnostic cut-off values and their performance for the evaluation of liver fibrosis severity by SWE™ - NASH - SuperSonic Imagine

\*\*\*\* Preliminary Internal R&D Study "Etude RND SSI hépatopathies chroniques 29072018" - SuperSonic Imagine

# CLINICAL CASE 4: PATIENT WITH KNOWN TYPE II DIABETES AND CIRRHOSIS

Parameters	Results	Pathology Scoring
Attenuation	0.64 dB/cm/MHz	Probable S2 *
SoS	1595 m/s	Probable S0 **
Stiffness	10.41 kPa	F3/F4***
Viscosity	2.17 Pa.s	Probable A3 / S2 ****



\* Fujiwara et al., The b-mode image-guided ultrasound attenuation parameter accurately detects hepatic steatosis in chronic liver disease, *Ultrasound Med. & Biol.*, Vol. 00, No. 00, pp. 1-10, 2018

\*\* Dioguardi Burgio et al., Ultrasonic Adaptive Sound Speed Estimation for the Diagnosis and Quantification of Hepatic Steatosis: A Pilot Study, *Ultraschall Med.* 2018 Nov 5.

\*\*\* MKG.EC.100 (Rev B) - Aixplorer Aixplorer Ultimate\_Guide assessment fibrosis steatosis - Proposed diagnostic cut-off values and their performance for the evaluation of liver fibrosis severity by SWE™ - NASH - SuperSonic Imagine

\*\*\*\* Preliminary Internal R&D Study "Etude RND SSI hépatopathies chroniques 29072018" - SuperSonic Imagine