



# VitaLogik<sup>®</sup> 4000 Series


A New Era in Computerized Monitoring



The VitaLogik 4000 Series Patient monitoring system is a compact monitor that suits a wide range of patient applications. The product is available in two versions: a non invasive monitor - VitaLogik 4000 and an invasive monitor - VitaLogik 4500. The VitaLogik 4000 series offer solutions to both intermediate and critical care monitoring. These systems are easily adapted to a diverse range of monitoring requirements including: Adult, Pediatric, Neonatal and patient transportation. Each device can be customized to meet the user's needs. As computerized monitors, the VitaLogik 4000 series, allow the users to take advantage of data connectivity, embedded data storage capabilities and high resolution display.




Big Numbers Display in the highest resolution available by a 12.1" Monitor. Allows easy data recognition from greater distance.



Has the ability to interface with, and display parameters from, other medical devices. Example data contains ventilator information



Capability for standard interface communication protocols included. Allows optional connectivity for additional data, including long term archived patient information, other network sources such as Laboratory, Radiology PACS and host based Hospital Systems. Includes a unique, automatic alert for lab results.



The VitaLogik 4000 series are the only 12.1" Display Monitors that allow up to 6 days Full Disclosure storage within the device.




Unique "All ECG Leads" display which is activated by a single command.

**Flexibility:** User defined "profiles," easily adaptable to a wide variety of patient applications and user specified requirements.  
Up to 6 hours battery time with all the vital signs available and therefore offering you a real portable / transport monitor.

**Automatic event collector:** When an alarm occurs, the monitor will automatically store the triggering event together with 20 seconds of full disclosure waveforms, trend graph and vital signs chart information that occurred at that specific time.

**The Clearest Images:** 12.1" monitor with 1024x768 resolution, provides you the sharpest and clearest images even when all 12 waveforms are displayed simultaneously.



**Clinical Added Value:** Simultaneous, multi-parameter display is capable of providing waveforms and numeric data for ECG (up to 12 leads, with ST and arrhythmia analysis), SpO2, Respiration, optional Invasive BP and Cardiac Output, along with Temperature and NIBP numeric data.

**Cost Efficient:** High-end capabilities, merged with low maintenance and service costs, minimize the total cost of ownership over the product's lifetime. This combination provides the best value in today's market and gives you the best overall purchase value.

**Ease of Use:** Has an intuitive operation and functionality, with an easy input method and straight forward navigation. Touch-screen display available (Optional).


*Monitor case components available in two different color versions!*

# VitaLogik 4000 Series Joins Mennen Medical Monitoring Solutions

Since 1963, Mennen Medical® has been committed to developing and marketing reliable, durable and top-quality products. Understanding the importance of accuracy in its field and implementing a zero-error policy, Mennen Medical® develops first-class equipment designed to suit the specific clinical needs and the dynamics of health care work environments.

From the Envoy™ - Mennen Medical's sophisticated, high-end modular bedside patient monitoring system, through the VitaLogik 5000 series (Invasive and Non Invasive) Multi-Parameters Monitors, to the Smart Central Stations and all the way to the Monitor's Web Portal that allows clinicians to access patient data stored in monitors from any location where an Internet connection is available, we provide a fully connected solution. The VitaLogik 4000 series is now joining the family.

Mennen Medical's monitoring solutions share a common user interface and an outstanding industrial design. It can be tailored to your needs, which include the easy transfer of patient data from VitaLogik 4000 series to the entire Mennen Medical patient monitoring products. Our people are here to make it happen for you.



[www.mennenmedical.com](http://www.mennenmedical.com)

**Headquarters:**  
Mennen Medical Ltd.  
4 Ha-Yarden Street,  
Yavne, P.O. Box 102,  
Rehovot 76100, Israel  
Phone: +972-8-9323333  
Fax: +972-8-9328510

**United States:**  
Mennen Medical Corp.  
101 Witmer Rd.,  
Suite # 100  
Horsham, PA 19044, US  
Phone: (215) 259 1020  
Fax: (215) 675 6212

**United Kingdom:**  
Mennen Medical UK  
Unit 18 Avant Business Center  
21 Denbigh Road  
Milton Keynes MK1 1DT, UK  
Phone: +44 1908 646070  
Fax: +44 1908 646030